

Exp 8: Functional Organic Materials and Molecular Modeling Preparation of 1',3',3'-Trimethylspiro(chromene-2,2'-indoline)

### Purpose

The goal of the experiment is to prepare 1',3',3'-trimethylspiro(chromene-2,2'-indoline and explore the thermochromism properties. The crystalline product will be analyzed with X-ray crystallography and NMR spectroscopy.



# Textbook reference

Brown, Iverson, Anslyn, & Foote, Section 19.5.

### Learning Objectives

Perform a synthesis.

Purify a product by recrystallization.

Use a melting point determination to analyze the product of a reaction.

Use IR spectroscopy to confirm the product of the reaction.

Analyze x-ray crystallographic data of the product.

# Equipment

- Reflux set up
- Büchner funnel with a vacuum filtration setup
- Rotary evaporator
- Steam bath
- Pyrex test tube

- Boiling stone
- Aluminum block heater
  - 1-dram vial with a
  - small hole in the cap
  - 42 MHz SpinSolve
- NMR and tubes

# Chemicals

- 1,3,3-Trimethyl-2methyleneindoline
- Salicylaldehyde
- Absolute ethanol
- Reagent grade
- acetone
- 90% Ethanol
- Diphenyl ether
- CDCl<sub>3</sub>

#### Part 1: Synthesis

- 1. Prepare a solution of 3.3 g of 1,3,3-trimethyl-2-methyleneindoline and 2.1 g of salicylaldehyde in 15 mL of absolute ethanol.
- 2. Heat the reaction mixture at reflux for 1 hour.
- 3. Filter the hot mixture with vacuum filtration and wash the filtered solid with a small amount of cold ethanol (5 mL).
- 4. Recrystallize the crude minor product (solid material on the filter paper) with acetone.
- 5. Obtain a mass of the minor product and determine the yield.
- 6. Perform a melting point determination of the minor product and record the value.
- 7. Evaporate the combined filtrate/alcohol wash using a rotary evaporator.
- 8. Recrystallize the resulting crude major product from 90% ethanol using a steam bath. Be careful not to melt the product.
- 9. Wash the crystals with cold absolute ethanol to reduce the amount of pigment in the sample.
- 10. Obtain a mass of the major product and determine the yield.
- 11. Perform a melting point determination of the major product and record the value.

#### Part 2: Thermochromism

- 12. Introduce a sample of the major product (about 250 mg) into a Pyrex test tube and add 5 mL of diphenyl ether. The test tube should be no more than one quarter full.
- 13. Add a boiling stone and gently boil the mixture by contacting the tube to a heated aluminum block on a hotplate.
- 14. Record your observations.

#### Part 3: Crystal growing

- 15. Large crystals suitable for X-ray crystallography (see Part 5) can be grown by dissolving 10–20 mg of your major product in 1 mL of acetone in a 1 dram vial with a small hole in the cap.
- 16. Leave the solution to evaporate in your locker.
- 17. Record your observations.

#### Part 4: Analysis using NMR spectroscopy

- 18. Dissolve 25 mg of your purified product in 0.5 mL CDCl<sub>3</sub> and transfer to an NMR tube.
- 19. Use the 42 MHz SpinSolve NMR in the classroom to obtain NMR spectra for your isolated product as well as the 1,3,3-trimethyl-2-methyleneindoline starting material.
- 20. Assign the structures of your product from these spectra.
- 21. Include the spectra with your report.

### Part 5: X-ray crystallography

- 22. Include crystal structures of major and minor products in your report. Crystallographic data can be obtained here: Xray Crystallographic Data.
- 23. Include a, b, c, and alpha, beta, gamma values for the unit cell of the structures (Some of these are not obvious, and a careful look into the Laue point group will explain what alpha, beta, and gamma should be for that group).

### Report

- 1. Write a full curved-arrow mechanism for this reaction, including the thermochromic process.
- 2. How does the NMR of the product and starting material differ? Could you use one (or more) peaks to track the progress of your reaction by NMR spectroscopy?
- 3. The first solid isolated in this reaction is actually a minor product which arises from the reaction of a second molecule of 1,3,3-trimethyl-2-methyleneindoline with the major product shown above. Draw the structure of this minor product.
- 4. Briefly explain what thermochromism is.
- 5. Based on the outcome of the thermochromism experiment, deduce which of the structures should be lower in energy, and explain your reasoning.
- 6. Draw a diagram of a simple x-ray experiment. Be sure to include the x-ray source, the crystal and the diffraction patterns obtained.
- 7. Thermochromic materials can be used in write-read-erase-rewrite media. Using the product of this experiment, discuss how you could design a recording medium based thereupon.