

Common Ion Effect



Objectives

- Determine solubility of KHT in a solution.
- Investigate effect of common ion on KHT solubility.
- Determine K_{sp} for KHT.



Solubility Equilibria

- **Solubility Rules** - **qualitative** prediction about solubility
- **Solubility Expression** - **quantitative** prediction about solubility or, how much will dissolve)
- The solid component is in equilibrium with the dissolved ions.



- **Solubility Product Constant, K_{sp}** – similar to other equilibrium expressions (K_{c} , K_{a} , etc.)

$$K_{\text{sp}} = [\text{Ba}^{2+}] [\text{SO}_4^{2-}]$$

- **Solubility** – quantity that dissolves to form a saturated solution.
 - Units are g/L or mol/L.




K_{sp} and Solubility

- K_{sp} remains constant with constant temperature.
- Solubility is affected by pH and concentration of other ions.

The Common-Ion Effect


- Solubility of a slightly soluble salt is *decreased* by presence of second solute with a common ion (**F⁻ ions**, in the example, below)



- Addition of NaF (**F⁻ ions**) to a solution of CaF_2
 - Shifts the equilibrium to the **← LEFT**, forming more CaF_2
 - Reducing its solubility; hence, precipitation)

Potassium Bitartrate (KHT)



- Today, we will study the solubility of *KHT*, a slightly soluble salt, a.k.a., *Cream of Tartar* (a by-product of winemaking ; used in baking powders).
- We will calculate K_{sp} in the presence of varying concentrations of a common ion.



$$K_{sp} = [\text{K}^+] [\text{HT}^-]$$

Potassium Bitartrate, cont.

- To solve for K_{sp} , we need $[K^+]$ and $[HT^-]$
- HT^- is a weak acid. We will titrate it with a strong base (NaOH).



At the equivalence point,

$$\text{moles } OH^- = \text{moles } HT^- \text{ and } \text{moles } HT^- = \text{moles KHT}$$


... which will lead you to solubility of KHT



Question: How does one know when the equivalence point is reached?

Answer: *Use an acid-base indicator!*

- Phenolphthalein
 - Use 2-3 drops only
 - Colorless in acidic solutions
 - Pink in solutions $\text{pH} \geq 8$



Faint pink color

The point where the indicator changes color is called the ***end point*** and is the approximation of the equivalence point.

Procedure : Step 1

- Nine (9) possible mixtures of KCl and NaCl
- See TA for your assigned mixture
 - Each pair will make 2 mixtures !!!
 - ALL mixtures must be chosen !!!
 - No more than 3 pairs of students can use the same mixture!


Mixture	0.1 M KCl (mL)	0.1 M NaCl (mL)
1	0	100
2	20	80
3	30	70
4	40	60
5	50	50
6	60	40
7	70	30
8	80	20
9	100	0

Procedure, cont.


- Perform vacuum filtration
- Use 2 burets (label both!)
 - Buret 1 – KHT solution
 - Use to measure 25.0 mL into an Erlenmeyer flask.
 - Record exact volume!
 - Buret 2 – NaOH solution
 - Use for titration to add the KHT to the Erlenmeyer flask which contains 2-3 drops of phenolphthalein indicator.
- Do three titrations for every mixture

Before leaving lab, see your TA and enter the following into the Excel spreadsheet on the lab computer:

- **Names**
- **Mixture #**
- **Volume KHT**
- **Volume NaOH**



Vacuum Filtration



Hazards and Waste

- Dispose all the solution in the Aqueous Waste container