


Equilibrium and Le Chatelier's Principle

Objectives

- To determine extinction coefficient for $[\text{FeSCN}]^{2+}$ and $[\text{CoSCN}]^{2+}$
- To calculate equilibrium concentrations
- To calculate the equilibrium constant, K_c
- Explore LeChatelier's Principle

Chemical Equilibrium

- The state where the concentrations of all reactants and products remain constant with time.
- At equilibrium, the rate of the forward reaction is equal to the rate of the reverse reaction.



Equilibrium constant K_c

$$K_c = \frac{[\text{C}]^c [\text{D}]^d}{[\text{A}]^a [\text{B}]^b}$$

Equilibrium does NOT depend on the amount of pure solids or pure liquids.

Absorbance and Beer's Law

- **Absorbance**
 - Amount of light absorbed by a particular species
- **Beer's Law**
 - Relates **absorbance to concentration**

$$A = \epsilon b c$$


Diagram illustrating the Beer's Law equation $A = \epsilon b c$ with labels for each term:

- A**: Absorbance (no units) (indicated by a blue arrow)
- ϵ : Extinction coefficient ($\text{L/mol}\cdot\text{cm}$) (indicated by a green arrow)
- b**: Path length = 1 cm (indicated by an orange arrow)
- c**: Concentration (M) (indicated by a dark blue arrow)


The plot of absorbance vs. concentration is linear and the slope is equal to ϵb .

Part A: Determining Extinction Coefficient and Equilibrium Concentrations


- 1st reaction: complex ion formation



- 2nd reaction: Complex ion formation



Part A: Determining ϵ and Equilibrium Concentrations (cont'd)




- Use the dilution equation to calculate $[\text{FeSCN}]^{2+}$ and $[\text{CoSCN}]^{2+}$.
 - (See tables in the procedure)

$$\text{Conc.}_{(\text{start})} \times \text{Volume}_{(\text{start})} = \text{Conc.}_{(\text{final})} \times \text{Volume}_{(\text{final})}$$


- Plot
 - Absorbance vs. $[\text{FeSCN}]^{2+}$
 - Absorbance vs. $[\text{CoSCN}]^{2+}$
- Find the Extinction coefficient ϵ for each reaction

LabQuest 3 Colorimeter

- Use 470 nm (**Blue**) LED for both the Iron and the Cobalt Reaction



- **Samples go in *cuvettes***

- These cost money and can be damaged so handle carefully
- Each group should only need one *Cuvette* that you will pick up at the back lab bench
- Use **ONE** *KimWipe* tissue to dry *cuvettes* and/or remove fingerprints
- Make sure the arrow  is pointing through the clear side of the cuvette and not the side with the ridges.
- Use the glass pipettes to prepare your samples.



Part B: Determining K_c

- Clean your pipet with DI water.
- Record absorbance for all 7 solutions:
 - 3 at room temperature and 4 at ice bath temperature
- Use ϵ from the previous graph to calculate $[\text{FeSCN}^{2+}]_{\text{eq}}$
- Use the Table & $M_1V_1=M_2V_2$ equation to find $[\text{Fe}^{3+}]_{\text{initial}}$ and $[\text{SCN}^-]_{\text{initial}}$
- Use ICE table to find $[\text{Fe}^{3+}]_{\text{eq}}$ and $[\text{SCN}^-]_{\text{eq}}$
- Plug values into the equilibrium expression to find K_c

$$K_c = \frac{[\text{C}]^c [\text{D}]^d}{[\text{A}]^a [\text{B}]^b}$$

ICE Table


Initial – **C**hange – **E**quilibrium

WRITE THIS DOWN NOW !


	Fe^{3+}	SCN^-	$[\text{FeSCN}]^{2+}$
Initial conc. (mol/L)	Determined using $C_1V_1 = C_2V_2$		0.00
Change (mol/L)	Related to the stoichiometry of the balanced equation		Difference between E quilibrium and I nitial values
Equilibrium (mol/L)	Difference between I nitial and C hange values		Calculated using ϵ and absorbance measurements



Part C: Le Chatelier's Principle



- A system will shift equilibrium position when disturbed to counteract the effect of the disturbance.
- Potential Disturbances
 - Change in component concentration
 - Change in temperature
 - Change in pressure
- Today's experiment will use solutions of cobalt(II) chloride. The shift in the **equilibrium** is accompanied by a change in color.
- The reaction of hydrated cobalt (II) ion and hydrochloric acid



Light pink

Purple/blue

Hazards and Waste

- Be very careful when using concentrated HCl. It will **VAPORIZE** so keep it in the fume hood so you do not inhale it.
- All waste should be placed into Aqueous Waste Container.